

C-DISC, konvex/TS

Produktmkm.: mit Titanspikes, konvexe Deckplatte, flache Grundplatte

Material: PEEK, Titan (Spikes)


Lieferinf.: steril



LFC Sp. z o.o.

C-DISC
 PEEK

Art.-Nr. REF	Winkelung (°)	Breite × Tiefe (mm)	Höhe post./ant. (mm)
869DP-00151304	0	15 × 13	4
869DP-00151305	0	15 × 13	5
869DP-00151306	0	15 × 13	6
869DP-00151307	0	15 × 13	7
869DP-00151308	0	15 × 13	8
881DP-07151304	7	15 × 13	4/6
881DP-07151305	7	15 × 13	5/7
881DP-07151306	7	15 × 13	6/8
881DP-07151307	7	15 × 13	7/9
881DP-07151308	7	15 × 13	8/10
869DP-00171504	0	17 × 15	4
869DP-00171505	0	17 × 15	5
869DP-00171506	0	17 × 15	6
869DP-00171507	0	17 × 15	7
869DP-00171508	0	17 × 15	8
881DP-07171504	7	17 × 15	4/6
881DP-07171505	7	17 × 15	5/7
881DP-07171506	7	17 × 15	6/8
881DP-07171507	7	17 × 15	7/9
881DP-07171508	7	17 × 15	8/10
875DP-00191504	0	19 × 15	4
875DP-00191505	0	19 × 15	5
875DP-00191506	0	19 × 15	6
875DP-00191507	0	19 × 15	7
881DP-07191504	7	19 × 15	4/6
881DP-07191505	7	19 × 15	5/7
881DP-07191506	7	19 × 15	6/8



C-DISC, flach/TS

Produktmkm.: mit Titanspikes, flache Deck- und Grundplatte

Material: PEEK, Titan (Spikes)

Lieferinf.: steril



LfC Sp. z o.o.



Art.-Nr. REF	Winkelung (°)	Breite × Tiefe (mm)	Höhe post./ant. (mm)
880DP-07151304	7	15 × 13	4/6
880DP-07151305	7	15 × 13	5/7
880DP-07151306	7	15 × 13	6/8
880DP-07151307	7	15 × 13	7/9
880DP-07151308	7	15 × 13	8/10
880DP-07171504	7	17 × 15	4/6
880DP-07171505	7	17 × 15	5/7
880DP-07171506	7	17 × 15	6/8
880DP-07171507	7	17 × 15	7/9
880DP-07171508	7	17 × 15	8/10

