

C-DISC, flach

Material: PEEK
Lieferinf.: steril



LCF Sp. z o.o.

Art.-Nr. REF	Winkelung (°)	Breite x Tiefe (mm)	Höhe post. (mm)	Höhe ant. (mm)	Auflagefl. je Seite (mm ²)
871DP-00151304	0	15 x 13	4	4	100 / 109
871DP-00151305	0	15 x 13	5	5	100 / 109
871DP-00151306	0	15 x 13	6	6	100 / 109
871DP-00151307	0	15 x 13	7	7	100 / 109
871DP-00151308	0	15 x 13	8	8	100 / 109
871DP-00151309	0	15 x 13	9	9	100 / 109
871DP-00171504	0	17 x 15	4	4	149 / 159
871DP-00171505	0	17 x 15	5	5	149 / 159
871DP-00171506	0	17 x 15	6	6	149 / 159
871DP-00171507	0	17 x 15	7	7	149 / 159
871DP-00171508	0	17 x 15	8	8	149 / 159
871DP-00171509	0	17 x 15	9	9	149 / 159
882DP-07151304	7	15 x 13	4	6	114 / 126
882DP-07151305	7	15 x 13	5	7	114 / 126
882DP-07151306	7	15 x 13	6	8	114 / 126
882DP-07151307	7	15 x 13	7	9	114 / 126
882DP-07151308	7	15 x 13	8	10	114 / 126
882DP-07151309	7	15 x 13	9	11	114 / 126
882DP-07171504	7	17 x 15	4	6	180 / 192
882DP-07171505	7	17 x 15	5	7	180 / 192
882DP-07171506	7	17 x 15	6	8	180 / 192
882DP-07171507	7	17 x 15	7	9	180 / 192
882DP-07171508	7	17 x 15	8	10	180 / 192
882DP-07171509	7	17 x 15	9	11	180 / 192



C-DISC, konvex

Produktmkm.: konvexe Deckplatte, flache Grundplatte

Material: PEEK

Lieferinf.: steril



LFC Sp. z o.o.



Art.-Nr. REF	Winkelung (°)	Breite x Tiefe (mm)	Höhe post. (mm)	Höhe ant. (mm)	Auflagefl. je Seite (mm ²)
872DP-00151304	0	15 x 13	4	4	109 / 109
872DP-00151305	0	15 x 13	5	5	109 / 109
872DP-00151306	0	15 x 13	6	6	109 / 109
872DP-00151307	0	15 x 13	7	7	109 / 109
872DP-00151308	0	15 x 13	8	8	109 / 109
872DP-00151309	0	15 x 13	9	9	109 / 109
872DP-00171504	0	17 x 15	4	4	168 / 191
872DP-00171505	0	17 x 15	5	5	168 / 191
872DP-00171506	0	17 x 15	6	6	168 / 191
872DP-00171507	0	17 x 15	7	7	168 / 191
872DP-00171508	0	17 x 15	8	8	168 / 191
872DP-00171509	0	17 x 15	9	9	168 / 191
883DP-07151304	7	15 x 13	4	6	100 / 126
883DP-07151305	7	15 x 13	5	7	100 / 126
883DP-07151306	7	15 x 13	6	8	100 / 126
883DP-07151307	7	15 x 13	7	9	100 / 126
883DP-07151308	7	15 x 13	8	10	100 / 126
883DP-07151309	7	15 x 13	9	11	100 / 126
883DP-07171504	7	17 x 15	4	6	154 / 192
883DP-07171505	7	17 x 15	5	7	154 / 192
883DP-07171506	7	17 x 15	6	8	154 / 192
883DP-07171507	7	17 x 15	7	9	154 / 192
883DP-07171508	7	17 x 15	8	10	154 / 192
883DP-07171509	7	17 x 15	9	11	154 / 192

