




LfC Sp. z o.o.

Cs-DISC, flach

Produktmkm.: 3D-TitanMatrix
Material: Ti6Al4V ELI
Lieferinf.: steril

Cs-DISC

Art.-Nr. REF	Winkelung (°)	Breite x Tiefe (mm)	Höhe post. (mm)	Höhe ant. (mm)	Auflagefl. je Seite (mm ²)
870DT-00151304	0	15 x 13	4	4	78
870DT-00151305	0	15 x 13	5	5	78
870DT-00151306	0	15 x 13	6	6	78
870DT-00151307	0	15 x 13	7	7	78
870DT-00151308	0	15 x 13	8	8	78
870DT-00151309	0	15 x 13	9	9	78
870DT-00151310	0	15 x 13	10	10	78
870DT-00171504	0	17 x 15	4	4	108
870DT-00171505	0	17 x 15	5	5	108
870DT-00171506	0	17 x 15	6	6	108
870DT-00171507	0	17 x 15	7	7	108
870DT-00171508	0	17 x 15	8	8	108
870DT-00171509	0	17 x 15	9	9	108
870DT-00171510	0	17 x 15	10	10	108
882DT-07151304	7	15 x 13	4	6	78
882DT-07151305	7	15 x 13	5	7	78
882DT-07151306	7	15 x 13	6	8	78
882DT-07151307	7	15 x 13	7	9	78
882DT-07151308	7	15 x 13	8	10	78
882DT-07151309	7	15 x 13	9	11	78
882DT-07151310	7	15 x 13	10	12	78
882DT-07171504	7	17 x 15	4	6	109
882DT-07171505	7	17 x 15	5	7	109
882DT-07171506	7	17 x 15	6	8	109
882DT-07171507	7	17 x 15	7	9	109
882DT-07171508	7	17 x 15	8	10	109
882DT-07171509	7	17 x 15	9	11	109
882DT-07171510	7	17 x 15	10	12	109





LF C Sp. z o.o.

Cs-DISC, konvex

Produktmkm.: 3D-TitanMatrix, konvexe Deckplatte, flache Grundplatte

Material: Ti6Al4V ELI

Lieferinf.: steril

Cs-DISC

Art.-Nr. REF	Winkelung (°)	Breite x Tiefe (mm)	Höhe post. (mm)	Höhe ant. (mm)	Auflagefl. je Seite (mm ²)
872DT-00151304	0	15 x 13	4	4	80 / 78
872DT-00151305	0	15 x 13	5	5	80 / 78
872DT-00151306	0	15 x 13	6	6	80 / 78
872DT-00151307	0	15 x 13	7	7	80 / 78
872DT-00151308	0	15 x 13	8	8	80 / 78
872DT-00151309	0	15 x 13	9	9	80 / 78
872DT-00151310	0	15 x 13	10	10	80 / 78
872DT-00171504	0	17 x 15	4	4	110 / 108
872DT-00171505	0	17 x 15	5	5	110 / 108
872DT-00171506	0	17 x 15	6	6	110 / 108
872DT-00171507	0	17 x 15	7	7	110 / 108
872DT-00171508	0	17 x 15	8	8	110 / 108
872DT-00171509	0	17 x 15	9	9	110 / 108
872DT-00171510	0	17 x 15	10	10	110 / 108
883DT-07151304	7	15 x 13	4	6	82 / 80
883DT-07151305	7	15 x 13	5	7	82 / 80
883DT-07151306	7	15 x 13	6	8	82 / 80
883DT-07151307	7	15 x 13	7	9	82 / 80
883DT-07151308	7	15 x 13	8	10	82 / 80
883DT-07151309	7	15 x 13	9	11	82 / 80
883DT-07151310	7	15 x 13	10	12	82 / 80
883DT-07171504	7	17 x 15	4	6	111 / 109
883DT-07171505	7	17 x 15	5	7	111 / 109
883DT-07171506	7	17 x 15	6	8	111 / 109
883DT-07171507	7	17 x 15	7	9	111 / 109
883DT-07171508	7	17 x 15	8	10	111 / 109
883DT-07171509	7	17 x 15	9	11	111 / 109
883DT-07171510	7	17 x 15	10	12	111 / 109

